THUASNE

Ligaflex® (2450)

- fr** Sangles d'immobilisation claviculaire
- en** Clavicular immobilisation straps
- de** Gurte zur Ruhigstellung des Schlüsselbeins
- nl** Riemen voor immobilisatie van het sleutelbeen
- it** Cinghie di immobilizzazione clavicolare
- es** Correas de inmovilización clavicolar
- pt** Correias de imobilização clavicolar
- da** Stropper til immobilisering af kraveben
- fi** Solisluiden immobilisaatioihnat
- sv** Immobiliseringsremmar för nyckelbenet
- el** Ιμάντες ακινητοποίησης κλείδας
- cs** Ortéza klíční kosti
- pl** Pasy unieruchamiające obojczyk
- lv** Atslēgas kaula imobilizācijas siksnas
- lt** Raktikaulių imobilizavimo diržai
- et** Rangluu immobiliseerimisrihmad
- sl** Trakovi za imobilizacijo ključnice
- sk** Klavikulárne popruhy na znehybnenie
- hu** Kulcscsont rögzítő pántok
- bg** Ремъци за имобилизация на ключицата
- ro** Curele de imobilizare claviculară
- ru** Ремни для имобилизации ключицы
- hr** Trake za imobilizaciju ključne kosti
- zh** 锁骨固定束帶
- ar** أحزمة تثبيت الترقوة

Ligaflex® Junior (2640)

Taille unique / Universal size

THUASNE
120, rue Marius AUFAN
92300 Levallois-Perret - France

© Thuasne - 6166301 (2021-10) © Studio Caterin

www.thuasne.com

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Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Figure of 8 straps allowing stabilisation of the clavicles for post-traumatic use.

Composition

Polyurethane foam - polyester - polyamide - cotton - polyethylene.

Properties/Mode of action

Two dorsal cross-straps fastened onto an anchoring pad. The simultaneous tightening of the straps helps stabilise the clavicles.

Adjustable and trimmable straps according to patient’s size.

Indications

Non-surgical clavicle fractures.

Contraindications

Do not use the product if the diagnosis has not been confirmed.

Do not use in the event of known allergy to any of the components. Do not apply the product in direct contact with broken skin.

Do not use in the event of complex or surgical shoulder fractures or rib cage fractures.

Do not use in case of severe vascular conditions of the upper limb (thrombosis, lymphedema, arterial compression...).

Precautions

Verify the product’s integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

Wear the device above thin clothing.

It is recommended that a healthcare professional supervises the first application.

Strictly comply with your healthcare professional’s prescription and recommendations for use.

For hygiene, security and performance reasons, do not re-use the product for another patient.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation.

Do not bring the product into direct contact with greasy substances (ointment, cream, etc.).

Do not wear the product when driving a vehicle.

The use of a medical device by a child, should be done under the supervision of an adult or a healthcare professional.

Undesirable side-effects

This device can cause skin reactions (redness, itching, burns, blisters, etc.) or wounds of various degrees of severity.

Possible risks of clavicle fracture complications.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

The help of another person is required to fit this device.

- After having cut the straps to the patient’s size (along the corresponding marking), place the self-fastening tabs at the end of the straps.
- Centre the dorsal anchoring pad and position the straps around the collar bones and under the armpits.
- Pass the straps through the buckles, adjust the tightening and fasten the straps via the self-fastening tabs.

Care

Product can be washed in accordance with the conditions shown on this notice and on the label.
2640: Hand wash.
2450: Machine washable at 40°C (delicate programme). If possible use a washing net. Close the self-fastening tabs before washing. Do not dry clean. Do not tumble-dry. Do not iron. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.).

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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Thuasne SAS – SIREN/RCS Nanterre 542 091 186