



**External  
breast  
prostheses**

**Bras for use  
with breast  
prostheses**

**SERENA**


**Made-to-measure  
sleeves or stockings  
for lymphoedema**



Sleeves or stockings for lymphoedema, made to measure for the patient by  
 Computer-Assisted Made-to-Measure Knitting (CAMMK) :  
 15-20 mmHg - 20-36 mmHg - > 36 mmHg.



Single sleeve



Sleeve + Mitten  
(separate or all-in-one)




### The 3 +

- + PRECISE**  
CAMMK ensures that the patient's measurements are respected.
- + COMFORTABLE**  
The softness of the microfibres and the finishes tailored to the condition promote patient compliance.
- + EFFECTIVE**  
The presence of polyamide makes it possible to control the elasticity and tailor it to the specific characteristics of the lymphoedema.




Shoulder-cover



Mitten



Women's stocking  
+ anti-slip, closed foot



Compression pants

### VARIOUS FINISHES

	UPPER LIMBS Sleeves 15-20 mmHg 20-36 mmHg > 36 mmHg	LOWER LIMBS Stockings and tights 15-20 mmHg 20-36 mmHg > 36 mmHg
<b>Distal end</b>	<ul style="list-style-type: none"> <li>• Single sleeve</li> <li>• Single sleeve + mitten (separate or all-in-one)</li> <li>• Single sleeve + mitten (separate or all-in-one) + thumb</li> </ul>	<ul style="list-style-type: none"> <li>• Open foot</li> <li>• Closed toe</li> </ul>
<b>Proximal end</b>	<ul style="list-style-type: none"> <li>• Single braid</li> <li>• Braid + anti-slip</li> <li>• Anti-slip for men or women</li> <li>• Shoulder cover</li> <li>• Shoulder strap</li> </ul>	<ul style="list-style-type: none"> <li>• Without anti-slip</li> <li>• Anti-slip for men or women</li> <li>• Compression pants</li> </ul>

- Composition : elastane, polyamide microfibres.
- Colour : beige.

Prescribe : THUASNE LYMPHOLOGY MADE-TO-MEASURE  
 sleeve or stocking



Right arm       Left arm


15-20 mmHg     20-36 mmHg     > 36 mmHg

Finish :

- Lace + Anti-slip
- Braid + Anti-slip
- Braid only
- Sleeve only
- Sleeve with mitten and thumb
- Sleeve with mitten without thumb
- Mitten only without thumb
- Mitten only with thumb
- Shoulder cover
- Shoulder attachment

circumferences

length measurements



**For hand, draw the end of the mitten and shade the zone covered by the orthosis.**

i      h

i = length 1-h (wrist line-end of mitten)      h = hand circumference at the end of the mitten

Measurements expressed in cm (e.g. 12.4 cm)

**IT IS ESSENTIAL TO COMPLETE ALL THE MEASUREMENT BOXES REQUIRED FOR THE DESIRED PRODUCT**

Prescriber's name and, if applicable, name of the medical institution : .....

.....

Patient's name : .....

.....

RETAILER'S ID

Customer N° :      Date :



## Silicone external breast prostheses

### • THUASNE know-how :

**SERENA** is a trademark of the Thuasne Group.











This family company specialised in the manufacture of medical products was created in 1847.

In order to meet the various requirements encountered in medicine, Thuasne has one strong strategic desire : to innovate and continuously improve the quality of its products and services.

**SERENA** external breast prostheses provide a supplement to orthopaedic ranges, with highly specialised technical products.

- When a woman has undergone a mastectomy, in the majority of cases, an external breast prosthesis remains the reference product to help her regain balance and comfort.
- Natural confidence means, above all, regaining her figure. For this, **SERENA** offers a wide range of models and sizes discreetly tailored to the shapes of your clients. This diversity enables a good fit, whatever type of surgery has been undergone.
- **SERENA** prostheses have been specially designed to provide maximum comfort. Thanks to its anatomical design, **SERENA** faithfully reproduces the breast shape and can be adjusted to fit each body shape without being noticed.
- In the blink of an eye, your client will regain her femininity and the confidence to resume all her social, professional and sporting activities (including swimming) without thinking twice.
- Prosthesis for external use only. • Similar consistency to natural tissues. • Excellent durability.

### • Types :

SERENA TRIANGULAR	SERENA HEART	SERENA ROUND	SERENA OVAL (SYMMETRICAL)	SERENA ASYMMETRICAL
				
				
Size : 1 to 10.	Size : 2 to 14.	Size : 1 to 11.	Size : 1 to 13.	Size : 1 to 13.
<b>Triangular shape</b> Very easy to adapt. Respects the natural position of the breast. Ref. 9201	<b>Anatomical shape</b> Perfectly tailored for a breast that is full in the upper part and rounded in the lower part. Particularly suitable for large sizes. Ref. 9341	<b>Traditional shape</b> Specifically tailored for evenly rounded breasts. Very suitable for small sizes. Ref. 9431	<b>Shape with lateral axillary extension</b> Versions : <b>round cavity (RC)</b> or <b>oval cavity (OC)</b> . The choice should be made on the basis of the scar. Ref. 9301 (RC) - Ref. 9311 (OC)	<b>Anatomical shape with lateral extension</b> Versions : <b>right - left</b> . To be recommended following a mastectomy with a pronounced axillary hollow. Ref. 9371 (R) - Ref. 9381 (L)

### • Sizes :

CUPS	BRAS SIZES	US/UK : <b>32</b>		US/UK : <b>34</b>		US/UK : <b>36</b>		US/UK : <b>38</b>		US/UK : <b>40</b>		US/UK : <b>42</b>		US/UK : <b>44</b>		US/UK : <b>46</b>		US/UK : <b>48</b>		US/UK : <b>50</b>	
		EUR : <b>70</b>	FR : <b>85</b>	EUR : <b>75</b>	FR : <b>90</b>	EUR : <b>80</b>	FR : <b>95</b>	EUR : <b>85</b>	FR : <b>100</b>	EUR : <b>90</b>	FR : <b>105</b>	EUR : <b>95</b>	FR : <b>110</b>	EUR : <b>100</b>	FR : <b>115</b>	EUR : <b>105</b>	FR : <b>120</b>	EUR : <b>110</b>	FR : <b>125</b>	EUR : <b>115</b>	FR : <b>130</b>
	Under chest measurement (in cm)	68 - 72	Serena Sizes	73 - 77	Serena Sizes	78 - 82	Serena Sizes	83 - 87	Serena Sizes	88 - 92	Serena Sizes	93 - 97	Serena Sizes	98 - 102	Serena Sizes	103 - 107	Serena Sizes	108 - 112	Serena Sizes	113 - 117	Serena Sizes
<b>A</b>	Chest measurement (in cm)	82 - 84	<b>2</b>	87 - 89	<b>3</b>	92 - 94	<b>4</b>	97 - 99	<b>5</b>	102 - 104	<b>6</b>	107 - 109	<b>7</b>	102 - 104	<b>8</b>	117 - 119	<b>9</b>	122 - 124	<b>10</b>	127 - 129	<b>11</b>
<b>B</b>		84 - 86	<b>3</b>	89 - 91	<b>4</b>	94 - 96	<b>5</b>	99 - 101	<b>6</b>	104 - 106	<b>7</b>	109 - 111	<b>8</b>	104 - 106	<b>9</b>	119 - 121	<b>10</b>	124 - 126	<b>11</b>	129 - 131	<b>12</b>
<b>C</b>		86 - 88	<b>4</b>	91 - 93	<b>5</b>	96 - 98	<b>6</b>	101 - 103	<b>7</b>	106 - 108	<b>8</b>	111 - 113	<b>9</b>	106 - 108	<b>10</b>	121 - 123	<b>11</b>	126 - 128	<b>12</b>	131 - 133	<b>13</b>
<b>D</b>		88 - 90	<b>5</b>	93 - 95	<b>6</b>	98 - 100	<b>7</b>	103 - 105	<b>8</b>	108 - 110	<b>9</b>	113 - 115	<b>10</b>	108 - 110	<b>11</b>	123 - 125	<b>12</b>	128 - 130	<b>13</b>	133 - 135	<b>14</b>

### • Indications :

	SERENA ASYMMETRICAL	SERENA OVAL (SYMMETRICAL)	SERENA TRIANGULAR	SERENA HEART	SERENA ROUND
Radical mastectomy	●●				
Modified mastectomy	●●	●●			
Mastectomy			●●	●●	●●

### • Packaging :

The SERENA breast prosthesis box also contains a small vanity case in which the prosthesis can be discreetly kept and carried around.



# Bras for use with breast prostheses




Elegant bras the shapes and fabrics of which have been technically designed for optimum comfort complete the breast prostheses range.

## Cotton flower



### Cotton bra




Colour : white

- Wide shoulder straps
- Feminine lace
- Fresh, soft feel of cotton
- Gentle wide elasticsation under the chest for greater comfort
- Internal lining of the cups for better support
- Fully covering model

## Peach skin



### Microfibre bra



Colour : white

- Wide shoulder straps
- Feminine lace
- Microfibre weave, fits softly like a second skin
- Wide basque under the chest for better support
- Inside cotton pocket
- Fully covering model

SERENA BRAS	Ref.
COTTON model - Right pocket	9800
COTTON model - Left pocket	9801
Microfibre model - Right pocket	9810
Microfibre model - Left pocket	9811

SERENA BRAS			
US/UK	EUR	F	CUPS
32	70	85	A - B - C
34	75	90	A - B - C - D
36	80	95	A - B - C - D
38	85	100	A - B - C - D
40	90	105	A - B - C - D
42	95	110	B - C - D

# SERENA

